



SLC30A7 Antibody, Biotin conjugated

Product Code	CSB-PA836732LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEW0
Immunogen	Recombinant Human Zinc transporter 7 protein (162-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 7 (ZnT-7) (Solute carrier family 30 member 7) (Znt-like transporter 2), SLC30A7, ZNT7 ZNTL2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC30A7