



ZDHHC1 Antibody, Biotin conjugated

Product Code	CSB-PA837848LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTX9
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC1 protein (283-389AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable palmitoyltransferase ZDHHC1 (EC 2.3.1.225) (DHHC domain-containing cysteine-rich protein 1) (Zinc finger DHHC domain-containing protein 1) (DHHC-1) (Zinc finger protein 377), ZDHHC1, C16orf1 ZNF377
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZDHHC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.