





TRIB3 Antibody

Product Code	CSB-PA847220LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RU7
Immunogen	Recombinant Human Tribbles homolog 3 protein (152-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tribbles homolog 3 (TRB-3) (Neuronal cell death-inducible putative kinase) (SINK) (p65-interacting inhibitor of NF-kappa-B), TRIB3, C20orf97 NIPK SKIP3 TRB3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRIB3

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA847220LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA847220LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.