



TRIB3 Antibody, HRP conjugated

Product Code	CSB-PA847220LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RU7
Immunogen	Recombinant Human Tribbles homolog 3 protein (152-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tribbles homolog 3 (TRB-3) (Neuronal cell death-inducible putative kinase) (SINK) (p65-interacting inhibitor of NF-kappa-B), TRIB3, C20orf97 NIPK SKIP3 TRB3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRIB3