



MPP5 Antibody, Biotin conjugated

Product Code	CSB-PA847626LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3R9
Immunogen	Recombinant Human MAGUK p55 subfamily member 5 protein (5-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MAGUK p55 subfamily member 5, MPP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.