




YTHDC1 Antibody

Product Code	CSB-PA853467LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MU7
Immunogen	Recombinant Human YTH domain-containing protein 1 protein (89-153AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	YTH domain-containing protein 1 (Splicing factor YT521) (YT521-B), YTHDC1, KIAA1966 YT521
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	YTHDC1

Image



Western Blot

Positive WB detected in: Mouse lung tissue, U251 whole cell lysate

All lanes: YTHDC1 antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 85, 83 kDa

Observed band size: 85 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.