



YTHDC1 Antibody, HRP conjugated

Product Code	CSB-PA853467LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MU7
Immunogen	Recombinant Human YTH domain-containing protein 1 protein (89-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	YTH domain-containing protein 1 (Splicing factor YT521) (YT521-B), YTHDC1, KIAA1966 YT521
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	YTHDC1