



RAB3IP Antibody, HRP conjugated

Product Code	CSB-PA853488LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QF0
Immunogen	Recombinant Human Rab-3A-interacting protein (1-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab-3A-interacting protein (Rab3A-interacting protein) (Rabin-3) (SSX2-interacting protein), RAB3IP, RABIN8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB3IP