




SLC29A1 Antibody

Product Code	CSB-PA857875LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99808
Immunogen	Recombinant Human Equilibrative nucleoside transporter 1 protein (228-291AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:2000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Equilibrative nucleoside transporter 1 (Equilibrative nitrobenzylmercaptapurine riboside-sensitive nucleoside transporter) (Equilibrative NBMPR-sensitive nucleoside transporter) (Nucleoside transporter, es-type) (Solute carrier family 29 member 1), SLC29A1, ENT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC29A1

Image



Western Blot

Positive WB detected in: Mouse skeletal muscle tissue


All lanes: KIR2DL3 antibody at 2.7µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 59 kDa

Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA857875LA01HU at dilution of 1:100



Immunofluorescent analysis of U251 cells using CSB-PA857875LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.