



SLC29A1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA857875LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99808 |
| Immunogen | Recombinant Human Equilibrative nucleoside transporter 1 protein (228-291AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Equilibrative nucleoside transporter 1 (Equilibrative nitrobenzylmercaptapurine riboside-sensitive nucleoside transporter) (Equilibrative NBMPR-sensitive nucleoside transporter) (Nucleoside transporter, es-type) (Solute carrier family 29 member 1), SLC29A1, ENT1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SLC29A1 |