



TNRC6C Antibody, HRP conjugated

Product Code	CSB-PA862070LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCJ0
Immunogen	Recombinant Human Trinucleotide repeat-containing gene 6C protein (157-421AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trinucleotide repeat-containing gene 6C protein, TNRC6C, KIAA1582
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TNRC6C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.