



PARD6G Antibody, HRP conjugated

Product Code	CSB-PA863155LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYG4
Immunogen	Recombinant Human Partitioning defective 6 homolog gamma protein (249-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Partitioning defective 6 homolog gamma (PAR-6 gamma) (PAR6D), PARD6G, PAR6G
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PARD6G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.