



MED17 Antibody, HRP conjugated

Product Code	CSB-PA865147LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVC6
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 17 protein (411-651AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 17 (Activator-recruited cofactor 77 kDa component) (ARC77) (Cofactor required for Sp1 transcriptional activation subunit 6) (CRSP complex subunit 6) (Mediator complex subunit 17) (Thyroid hormone receptor-associated protein complex 80 kDa component) (Trap80) (Transcriptional coactivator CRSP77) (Vitamin D3 receptor-interacting protein complex 80 kDa component) (DRIP80), MED17, ARC77 CRSP6 DRIP77 DRIP80 TRAP80
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MED17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.