



ING3 Antibody, Biotin conjugated

Product Code	CSB-PA865171LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXR8
Immunogen	Recombinant Human Inhibitor of growth protein 3 protein (51-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of growth protein 3 (p47ING3), ING3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ING3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.