




IQCG Antibody

Product Code	CSB-PA867116LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H095
Immunogen	Recombinant Human IQ domain-containing protein G protein (1-134AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IQ domain-containing protein G, IQCG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IQCG

Image



Western Blot

Positive WB detected in: Mouse heart tissue,
Mouse skeletal muscle tissue


All lanes: IQCG antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52, 41 kDa

Observed band size: 52, 4 kDa



Immunofluorescent analysis of HepG2 cells
using CSB-PA867116LA01HU at dilution of
1:100 and Alexa Fluor 488-conjugated
AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.