





PPP1R12C Antibody

Product Code	CSB-PA871627LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZL4
Immunogen	Recombinant Human Protein phosphatase 1 regulatory subunit 12C protein (317-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1 regulatory subunit 12C (Protein phosphatase 1 myosin-binding subunit of 85 kDa) (Protein phosphatase 1 myosin-binding subunit p85), PPP1R12C, LENG3 MBS85
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1R12C

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA871627LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA871627LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)