



PPP1R12C Antibody, Biotin conjugated

Product Code	CSB-PA871627LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZL4
Immunogen	Recombinant Human Protein phosphatase 1 regulatory subunit 12C protein (317-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1 regulatory subunit 12C (Protein phosphatase 1 myosin-binding subunit of 85 kDa) (Protein phosphatase 1 myosin-binding subunit p85), PPP1R12C, LENG3 MBS85
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1R12C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.