



SLC52A3 Antibody, Biotin conjugated

Product Code	CSB-PA873618LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ40
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 3 protein (159-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 3 (Riboflavin transporter 2) (hRFT2), SLC52A3, C20orf54 RFT2 RFVT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC52A3