



SLC52A3 Antibody, Biotin conjugated

Product Code	CSB-PA873618LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ40
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 3 protein (159-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 3 (Riboflavin transporter 2) (hRFT2), SLC52A3, C20orf54 RFT2 RFVT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC52A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.