



ZC4H2 Antibody, HRP conjugated

Product Code	CSB-PA873628LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQZ6
Immunogen	Recombinant Human Zinc finger C4H2 domain-containing protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger C4H2 domain-containing protein (Hepatocellular carcinoma-associated antigen 127), ZC4H2, HCA127 KIAA1166
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC4H2