



ZC4H2 Antibody, Biotin conjugated

Product Code	CSB-PA873628LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQZ6
Immunogen	Recombinant Human Zinc finger C4H2 domain-containing protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger C4H2 domain-containing protein (Hepatocellular carcinoma-associated antigen 127), ZC4H2, HCA127 KIAA1166
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC4H2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.