



SMURF1 Antibody, Biotin conjugated

Product Code	CSB-PA875720LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCE7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase SMURF1 protein (198-374AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase SMURF1 (hSMURF1) (EC 2.3.2.26) (HECT-type E3 ubiquitin transferase SMURF1) (SMAD ubiquitination regulatory factor 1) (SMAD-specific E3 ubiquitin-protein ligase 1), SMURF1, KIAA1625
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SMURF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.