



PGAP2 Antibody, Biotin conjugated

Product Code	CSB-PA883397LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHH9
Immunogen	Recombinant Human Post-GPI attachment to proteins factor 2 protein (45-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Post-GPI attachment to proteins factor 2 (FGF receptor-activating protein 1), PGAP2, FRAG1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PGAP2