



RACGAP1 Antibody, HRP conjugated

Product Code	CSB-PA884432LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0H5
Immunogen	Recombinant Human Rac GTPase-activating protein 1 protein (86-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rac GTPase-activating protein 1 (Male germ cell RacGap) (MgcRacGAP) (Protein CYK4 homolog) (CYK4) (HsCYK-4), RACGAP1, KIAA1478 MGCRACGAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RACGAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.