



SLC28A3 Antibody, HRP conjugated

Product Code	CSB-PA884504LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAS3
Immunogen	Recombinant Human Solute carrier family 28 member 3 protein (1-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 28 member 3 (Concentrative Na(+)-nucleoside cotransporter 3) (CNT 3) (hCNT3), SLC28A3, CNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC28A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.