



SLC28A3 Antibody, FITC conjugated

Product Code	CSB-PA884504LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAS3
Immunogen	Recombinant Human Solute carrier family 28 member 3 protein (1-102AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 28 member 3 (Concentrative Na(+)-nucleoside cotransporter 3) (CNT 3) (hCNT3), SLC28A3, CNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC28A3