




METTL14 Antibody

Product Code	CSB-PA884519LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCE5
Immunogen	Recombinant Human N6-adenosine-methyltransferase subunit METTL14 protein (24-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	METTL14

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate

All lanes: METTL14 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.