



METTL14 Antibody, HRP conjugated

Product Code	CSB-PA884519LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCE5
Immunogen	Recombinant Human N6-adenosine-methyltransferase subunit METTL14 protein (24-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N6-adenosine-methyltransferase non-catalytic subunit (Methyltransferase-like protein 14) (hMETTL14), METTL14, KIAA1627
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	METTL14