




ALKBH4 Antibody

Product Code	CSB-PA885759LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXW9
Immunogen	Recombinant Human Alpha-ketoglutarate-dependent dioxygenase alkB homolog 4 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-ketoglutarate-dependent dioxygenase alkB homolog 4 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 4), ALKBH4, ABH4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ALKBH4

Image



Western Blot

Positive WB detected in: Mouse spleen tissue,

Mouse lung tissue, Mouse kidney tissue


All lanes: ALKBH4 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 34, 5 kDa

Observed band size: 34 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA885759LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.