



GGA1 Antibody, HRP conjugated

Product Code	CSB-PA890911LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJY5
Immunogen	Recombinant Human ADP-ribosylation factor-binding protein GGA1 protein (1-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-binding protein GGA1 (Gamma-adaptin-related protein 1) (Golgi-localized, gamma ear-containing, ARF-binding protein 1), GGA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GGA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.