



PDZRN3 Antibody, Biotin conjugated

Product Code	CSB-PA891797LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPQ7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase PDZRN3 protein (133-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase PDZRN3 (EC 2.3.2.27) (Ligand of Numb protein X 3) (PDZ domain-containing RING finger protein 3) (RING-type E3 ubiquitin transferase PDZRN3) (Semaphorin cytoplasmic domain-associated protein 3) (Protein SEMACAP3), PDZRN3, KIAA1095 LNX3 SEMCAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PDZRN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.