



NOTCH3 Antibody, Biotin conjugated

Product Code	CSB-PA892153LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM47
Immunogen	Recombinant Human Neurogenic locus notch homolog protein 3 protein (1301-1496AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain], NOTCH3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOTCH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.