



TACC3 Antibody, HRP conjugated

Product Code	CSB-PA896927LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6A5
Immunogen	Recombinant Human Transforming acidic coiled-coil-containing protein 3 protein (133-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming acidic coiled-coil-containing protein 3 (ERIC-1), TACC3, ERIC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TACC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.