




RTCB Antibody

Product Code	CSB-PA897546LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3I0
Immunogen	Recombinant Human tRNA-splicing ligase RtcB homolog protein (52-229AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA-splicing ligase RtcB homolog (EC 6.5.1.3), RTCB, C22orf28
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RTCB

Image



Western Blot

Positive WB detected in: NIH/3T3 whole cell lysate, HEK293 whole cell lysate


All lanes: RTCB antibody at 3μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56 kDa


Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA897546LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA897546LA01HU at dilution of 1:100



Immunofluorescent analysis of Hela cells using CSB-PA897546LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)