



TRAF6 Antibody, HRP conjugated

Product Code	CSB-PA024154LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4K3
Immunogen	Recombinant Human TNF receptor-associated factor 6 protein (1-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TNF receptor-associated factor 6 (EC 2.3.2.27) (E3 ubiquitin-protein ligase TRAF6) (Interleukin-1 signal transducer) (RING finger protein 85) (RING-type E3 ubiquitin transferase TRAF6), TRAF6, RNF85
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRAF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.