



L1 Antibody, Biotin conjugated

Product Code	CSB-PA08349D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03101
Immunogen	Recombinant Human papillomavirus type 16 Major capsid protein L1 protein (6-471AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 16
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein L1, L1
Immunogen Species	Human papillomavirus type 16
Research Area	Microbiology
Target Names	L1