



L1 Antibody, Biotin conjugated

Product Code	CSB-PA08349D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03101
Immunogen	Recombinant Human papillomavirus type 16 Major capsid protein L1 protein (6-471AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 16
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein L1, L1
Immunogen Species	Human papillomavirus type 16
Research Area	Microbiology
Target Names	L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.