



Recombinant Human metapneumovirus Fusion glycoprotein F0 (F), partial

Product Code	CSB-MP751041HDAM1
Abbreviation	F
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6WB98
Product Type	Recombinant Protein
Immunogen Species	Human metapneumovirus (strain CAN97-83) (HMPV)
Purity	>85% (SDS-PAGE)
Sequence	KESYLEESCSTITEGYLSVLRGTWYTNVFTLEVGDVENLTCSDGPSLIKTELDL TKSALR ELKTVSADQLAREEQIENPRQSR
Source	Mammalian cell
Target Names	F
Protein Names	Recommended name: Fusion glycoprotein F0 Short name= Protein FCleaved into the following 2 chains: 1. Fusion glycoprotein F2 2. Fusion glycoprotein F1
Expression Region	20-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.