



Recombinant Human HLA class I histocompatibility antigen, alpha chain G (HLA-G), partial

Product Code	CSB-EP010509HU1-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17693
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>GSHSMRYFSAAVSRPGRGEPRIAMGYVDDTQFVRFSDSACPRMEPRAPW VEQEGPEYW EEETRNTKAHAQTDRMNLQTLRGYYNQSEASSHTLQWMIGCDLGSDGRLLR GYEQYAYDG KDYLALNEDLRSWTAADTAAQISKRKCEANVAEQRRAYLEGTCVEWLHRYL ENKEMLQ RADPPKTHVTHHPVFDYEATLRCWALGFYPAEIILTWQRDGEDQTQDVELVET RPAGDGT FQKWA AVV VPSGEEQRYTCHVQHEGLPEPLMLRWKQSSSLPTIPI</p>
Source	E.coli
Target Names	HLA-G
Protein Names	<p>Recommended name: HLA class I histocompatibility antigen, alpha chain G Alternative name(s): HLA G antigen MHC class I antigen G</p>
Expression Region	25-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>