



# Recombinant Human Taste receptor type 2 member 46 (TAS2R46), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP023163HU1   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P59540  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | WSFESLENKPVFMFCEAIAFSYPSTHPFILIWGNKCLKQTFLSVLWHVRYWVK<br>GEKPSSS  |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | TAS2R46   |
| <b>Protein Names</b>     | Recommended name: Taste receptor type 2 member 46 Short name=<br>T2R46Alternative name(s): Taste receptor type 2 member 54 Short name=<br>T2R54   |
| <b>Expression Region</b> | 250-309aa   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |