



GAPDH Antibody, Biotin conjugated

Product Code	CSB-PA00029D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04406
Immunogen	Recombinant Human Glyceraldehyde-3-phosphate dehydrogenase protein (2-335AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyceraldehyde-3-phosphate dehydrogenase (GAPDH) (EC 1.2.1.12) (Peptidyl-cysteine S-nitrosylase GAPDH) (EC 2.6.99.-), GAPDH, GAPD
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GAPDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.