




ABCG1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001080LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P45844 |
| Immunogen | Recombinant Human ATP-binding cassette sub-family G member 1 protein (84-222AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | ATP-binding cassette sub-family G member 1 (ATP-binding cassette transporter 8) (White protein homolog), ABCG1, ABC8 WHT1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | ABCG1 |

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA001080LA01HU at dilution of 1:100