




ABCG1 Antibody

Product Code	CSB-PA001080LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45844
Immunogen	Recombinant Human ATP-binding cassette sub-family G member 1 protein (84-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family G member 1 (ATP-binding cassette transporter 8) (White protein homolog), ABCG1, ABC8 WHT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ABCG1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA001080LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.