



ABCG1 Antibody, FITC conjugated

Product Code	CSB-PA001080LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45844
Immunogen	Recombinant Human ATP-binding cassette sub-family G member 1 protein (84-222AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family G member 1 (ATP-binding cassette transporter 8) (White protein homolog), ABCG1, ABC8 WHT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ABCG1