



ADRB1 Antibody, FITC conjugated

Product Code	CSB-PA001391LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08588
Immunogen	Recombinant Human Beta-1 adrenergic receptor protein (401-477AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1 adrenergic receptor (Beta-1 adrenoreceptor) (Beta-1 adrenoceptor), ADRB1, ADRB1R B1AR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ADRB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.