



ALDH1A3 Antibody, FITC conjugated

Product Code	CSB-PA001567LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47895
Immunogen	Recombinant Human Aldehyde dehydrogenase family 1 member A3 protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aldehyde dehydrogenase family 1 member A3 (EC 1.2.1.5) (Aldehyde dehydrogenase 6) (Retinaldehyde dehydrogenase 3) (RALDH-3) (RaLDH3), ALDH1A3, ALDH6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ALDH1A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.