




APTX Antibody

Product Code	CSB-PA001955LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z2E3
Immunogen	Recombinant Human Aprataxin protein (92-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aprataxin (EC 3.1.11.7) (EC 3.1.12.2) (Forkhead-associated domain histidine triad-like protein) (FHA-HIT), APTX, AXA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	APTX

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA001955LA01HU at dilution of 1:100