



ARHGDIB Antibody, Biotin conjugated

Product Code	CSB-PA002038LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52566
Immunogen	Recombinant Human Rho GDP-dissociation inhibitor 2 protein (2-86AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GDP-dissociation inhibitor 2 (Rho GDI 2) (Ly-GDI) (Rho-GDI beta), ARHGDIB, GDIA2 GDID4 RAP1GN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGDIB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.