



# ARHGDIB Antibody, Biotin conjugated

Product Code	CSB-PA002038LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52566
Immunogen	Recombinant Human Rho GDP-dissociation inhibitor 2 protein (2-86AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GDP-dissociation inhibitor 2 (Rho GDI 2) (Ly-GDI) (Rho-GDI beta), ARHGDIB, GDIA2 GDID4 RAP1GN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGDIB