




ATXN1L Antibody

Product Code	CSB-PA002440LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7T5
Immunogen	Recombinant Human Ataxin-1-like protein (183-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ataxin-1-like (Brother of ataxin-1) (Brother of ATXN1), ATXN1L, BOAT BOAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ATXN1L

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA002440LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.